

Meril

More to Life

The Next Era of Urology

merillife.com



Cardiovascular

Orthopedics

Diagnostics

Endo-Surgery

Robotics

More to Life

Meril is a global medical device company dedicated to the design and development of novel, clinically relevant, state-of-the-art and best-in-class devices to alleviate human suffering and improve the quality of life. We span a broad operational canvas ranging from vascular interventional devices to orthopedics, in-vitro diagnostics and endo-surgery.

We share an enduring perseverance to advance healthcare solutions, for more patients to live longer and healthier lives. We have a strong commitment towards R&D and adhere to the best quality standards in manufacturing, scientific communication and distribution known today.

Established in 2006, Meril's launch was in line with the healthcare diversification plan of Bilakhia Group. Located in Vapi, 150km north of Mumbai, and ensconced within the serene hill-sides of Chala town, Meril sits proud in its ultra modern manufacturing facilities.

Where research becomes the basis for discovering new ways of alleviating human suffering and empowering an individual with good health.

 **Research
Empowering
Individual**





Scope Defined with Clarity

Urology covers a broad spectrum of conditions affecting the urinary system and male reproductive health. From kidney stones to prostate care, the scope of urology extends to diagnosing, managing, and treating diverse urological disorders with cutting-edge solutions.



Access the Gateway to Treatment

Timely access to urological care is crucial for better patient outcomes. With innovations in minimally invasive techniques and robotic-assisted procedures, urology ensures patients receive precise, effective, and less invasive treatments faster than ever.





Diagnose Early & Accurately

Accurate diagnosis is the foundation of effective treatment in urology. With state-of-the-art imaging, diagnostics, and AI-driven assessments, specialists can detect conditions at an early stage, ensuring personalized and proactive patient care.



Treat Smarter, Heal Faster

Accurate diagnosis is the foundation of effective treatment in urology. With state-of-the-art imaging, diagnostics, and AI-driven assessments, specialists can detect conditions at an early stage, ensuring personalized and proactive patient care.

01 Consumable

Catheter

- 11 Deflure™ Ureteral Balloon Dilation Catheter
- 13 SupraSure™ Suprapubic Cystostomy Catheter
- 15 Rescath™ General-Purpose Ureteral Catheter

Stone Basket

- 19 Ri-Catch™ Urinary Stone Retrieval Basket {F}
- 21 Ri-Catch™ + Urinary Stone Retrieval Basket {R}
- 23 Messtone™ Urological Retrieval Coil

Sheath

- 27 Provure™ Ureteral Access Sheath{R}
- 29 Provure P™ + Ureteral Access Sheath{Y}
- 31 Provure Pro™ Ureteral Access Sheath{YB}
- 33 Glizer™ Amplatz Sheath

Stent

- 37 UretrixPro™ Ureteral Stent Package
- 39 Uretrix™ Ureteral Stent
- 41 Uretrix™ + Hydrophilic Ureteral Stent
- 43 ReMosten™ Mono Stent

Guidewire

- 47 Zebcor™ Zebra guidewire
- 49 Ribrid™ Hybrid Guidewire
- 51 Hydroease™ Hydrophilic Guidewire
- 53 Meriskate™ PTFE Guidewire

Dilator

- 57 Dilato™ Fascial Dilator Set
- 59 Dilatrix™ Urethral Dilator
- 61 EasyDi™ Ureteral Dilator

02 Scope

65 ri-Uniview™ Flexible Fibreoptic Ureteroscope

67 riUroSnap™ Ureteroscope System

03 Laser

71 MeriFRAG+™ Holmium Laser Portable System

04 Special Category

73 ProNeph™ Percutaneous Nephrolithotomy (PCNL) Package

The Next Era of Urology



Catheter

Precision Flow, Seamless Motion

Our range of catheters ensures optimal flow and seamless performance in urological care. Designed for efficiency, comfort, and reliability, our solutions empower both patients and healthcare professionals to move forward with confidence.

DeFlure™
Ureteral Balloon Dilatation Catheter

SupraSure™
Suprapubic Cystostomy Catheter

Rescath™
General-Purpose Ureteral Catheter





DeflureTM

Ureteral Balloon Dilatation Catheter

Features

- Dual markers for precise balloon positioning Reduces the need for prior dilation or guidewire manipulation.
- Ensures rapid inflation/deflation, crucial in time-sensitive urologic procedures.
- Pebax Material for optimal flexibility and pushability, smooth navigation and precise dilation.

Size Matrix

Catheter	Balloon Length	Catheter Length	Inflated Balloon Diameter
4 FR	60 mm	850 mm	5 mm
4 FR	60 mm	650 mm	

Use

The catheter dilates ureteral strictures or prepares the ureter for procedures.



SupraSureTM

Suprapubic Cystostomy Catheter

Features

- Chrysanthemum tip design can effectively prevent slip
- Handle design for easy grip
- Silicone remains soft at body temperature, ensuring enhanced patient comfort during extended use.

Size Matrix

Catheter size	18 Fr Length: 23 Cm
---------------	---------------------

Use

It provides bladder drainage via percutaneous Malecot catheter placement.



RescathTM

General-Purpose Ureteral Catheter

Features


- Crafted from soft PVC for ideal flexibility and lasting strength.
- Enables smooth navigation and precise contrast delivery in urology procedures.
- Tapered tip reduces ureteral trauma and improves patient comfort.

Size Matrix

5 FR Length	6 FR Length
70 CM	70 CM

Use

It's a flexible tube for insertion into the ureters via cystoscope or nephroscope.



The Next Era of **Urology**

Stone Basket

Secure Grip, Precise Extraction

Just like an athlete's perfect catch, our stone baskets ensure a firm yet flexible grip for precise kidney stone retrieval. With shape memory and flexible properties, they minimize mucosal trauma, bleeding, and injuries while maximizing procedural efficiency—empowering surgeons with confidence for every extraction.

Ri-Catch[™]

Urinary Stone Retrieval Basket (Flexible)

Ri-Catch +[™]

Urinary Stone Retrieval Basket (Rigid)

Messtone[™]

Urological Retrieval Coil





Ri-CatchTM

Urinary Stone Retrieval Basket (Flexible)

Features

- Shape-memory design quickly restores basket form after sharp bends.
- Tipless basket reduces trauma and enhances clear visualization.
- Nitinol wire offers strong shape memory and twist resistance for basket stability.
- PI core boosts rotation and control for precise stone retrieval.


Size Matrix

Basket Size(Flexible)

L:90CM	L:110CM	L:120CM
(3 F) (4 wire)	1.7FR (4 wire) ID – 0.38 mm, OD – 0.56 mm	(1.7) (4 wire)
		(2.2 FR) (4 wire)
(1.9FR) (3 wire)	2.2FR (4 wire) ID – 0.45 mm, OD – 0.73 mm	(3 F) (4 wire)
(2.2FR) (3 wire)		(1.9FR) (3 wire)
		(2.2FR) (3 wire)

Use

It removes urinary stones during endoscopic procedures.



Ri-Catch +TM

Urinary Stone Retrieval Basket (Rigid)

Features

- Flexible shape-memory basket quickly returns to form after bends.
- Tipless design minimizes trauma and ensures clear visualization.
- It is used with rigid endoscope

Size Matrix

Basket Size(Rigid): 4.5 FR (4 wire), length: 55 cm

Use

It removes urinary stones during rigid endoscopic procedures.



MesstoneTM

Urological Retrieval Coil

Features

- It allows secure extraction of stone pieces.
- Can resist against unwanted holmium laser light energy.
- Facilitates fast exchanges of scopes with quick release handle.
- Enhances stone fragmentation, boosts stone-free rate, and reduces need for additional therapy.

Size Matrix

Stone coil: 3FR x Length: 115 CM

Use

It is used endoscopically to capture and remove urinary tract stones and objects.



The Next Era of **Urology**

Sheath

Seamless Access, Superior Strength

Just like a javelin's aerodynamic design ensures precision and control, our sheaths provide kink resistance and pressure resistance for uncompromised access and stability. Engineered for effortless insertion and smooth navigation, they minimize trauma while maximizing procedural efficiency, empowering surgeons with precision and confidence.

ProvureTM
Ureteral Access Sheath

Provure ProTM
Ureteral Access Sheath - YB Bendable

ProvureTM +
Ureteral Access Sheath - Y Type

GlizerTM
Amplatz Sheath



ProvureTM

Ureteral Access Sheath

Features


- Hydrophilic coating smooths the sheath upon contact with urine.
- Enables multiple scope entries without re-dilation, reducing OR time and ensuring smooth tool operation.
- Maintains lumen integrity during manipulation.

Size Matrix

L:26 CM	9/11 FR		0/12 FR
L:35 CM	2/14 FR	9.5/11.5 FR	10.5/12.5 FR
L:45 CM	9.5/11.5 FR	10.5/12.5 FR	12/14 FR

Use

It helps as a working channels for the instruments during endoscopic urinary procedures.



ProvureTM +

Ureteral Access Sheath - Y Type

Features


- Facilitates continuous irrigation and efficient fluid management.
- Enables simultaneous scope access and suction.
- Improves visibility by rapidly clearing stone fragments and debris
- Y connector available for negative suction

Size Matrix

L:26 CM	L:40 CM
9/11 FR	14/16FR
10/12 FR	

Use

It facilitates instrument passage and Suction in endoscopic urinary procedures.



ProvureTM Pro

Ureteral Access Sheath - YB Bendable

Features

- Enhances access to lower pole and complex kidney calyces
- Maintains lumen integrity while bending—ensuring consistent scope passage
- Y connector available for negative suction with bendable tip

Size Matrix

L:26 CM	9/11 FR		10/12 FR
L:35 CM	9.5/11.5 FR	10.5/12.5 FR	12/14 FR
L:45 CM	9.5/11.5 FR	10.5/12.5 FR	12/14 FR
L:46 CM	13/15FR		

Use

It aids endoscope passage and fluid injection in urinary procedures.



GlizerTM

Amplatz Sheath

Features

- Tapered distal end allows smooth, atraumatic advancement over dilators or balloons.
- Amplatz sheath's kink-resistant PTFE boosts durability and reliability.

Size Matrix

L:15 CM	L:17 CM	L:20 CM	L:30 CM
18FR	18FR	12FR	18FR
		18FR	
20FR	20FR	24FR	20FR
		30 FR	
22FR	22FR	14FR	22FR
		20FR	
24FR	24FR	26FR	24FR
		16FR	
26FR	26FR	22FR	26FR
		28FR	
28FR	28FR		28FR

Use

It dilates tissue around a dilator to maintain a nephrostomy tract for instrument access.



Stents

Strength in Support, Precision in Placement

Just as an archer relies on skill and precision to hit the mark, surgeons need accuracy and control for successful interventions. Our stents are designed for effortless deployment, optimal positioning, and long-term reliability—empowering surgeons with confidence and ensuring seamless outcomes for patients.

Uretrix[™] Pro
Ureteral Stent Package

Uretrix[™] +
Hydrophilic Ureteral Stent

RiMosten[™]
Mono Stent

Uretrix[™]
Ureteral Stent



UretrixTM Pro

Ureteral Stent Package

Features

- Coating: With or without hydrophilic coating.
- In-dwelling Time: Less than 3 months.
- The ureteral stent is temperature-controlled, causing it to soften once inserted.
- Multiple holes with large lumens lead to improved drainage.
- Tapered tip is designed for easy advancement.

Size (Fr)

L:26 CM	L:28 CM	Guidewire
4.8 FR	4.8 FR	0.035 inch
5 FR	5 FR	
6 FR	6 FR	

Use

The ureteral stent provides drainage to relieve obstructions.



UretrixTM +

Hydrophilic Ureteral Stent

Features

- Reduces surface friction when wet, ensuring smooth insertion and removal.
- In-dwelling Time: Less than 3 months.
- Inner lumen design allows uninterrupted urine drainage, reducing risk of blockage.
- Multiple holes with large lumens lead to improved drainage.
- Low-friction surface reduces irritation and post-operative discomfort.

Size Matrix

Stent Size 26 CM	Stent Size 28 CM
4.8 FR	4.8 FR
5 FR	5 FR
6 FR	6 FR

Use

Hydrophilic ureteral stent provides temporary drainage from the ureteropelvic junction to the bladder.



UretrixTM

Ureteral Stent

Features

- Coating: With or without hydrophilic coating
- In-dwelling Time: Less than 3 months.
- The ureteral stent is temperature-controlled, causing it to soften once inserted.
- Multiple holes with large lumens lead to improved drainage.
- Tapered tip is designed for easy advancement.

Size Matrix

L:26 CM	L:28 CM
4.8 FR	4.8 FR
5 FR	5 FR
6 FR	6 FR

Use

Ureteral stent provides temporary drainage from the ureteropelvic junction to the bladder.



RiMostenTM

Mono Stent

Features

- Made from kink-resistant, durable polyurethane for reliable flexibility.
- In-dwelling Time: Less than 3 months.
- Holes in shaft lead to improved drainage
- Tapered tip is designed for easy advancement
- Ensures secure renal anchorage with minimal
- Ureteral and bladder irritation.

Size Matrix

10 FR x Length: 50 CM

Use

Mono stent is used for drainage and irrigation, with a pig tail tip to retain it in the renal pelvis.



Guidewire

Smooth Navigation, Absolute Control

Just like a skier glides smoothly through twists and turns, our guidewires ensure seamless movement through complex pathways. Featuring a super-flexible nitinol core that inhibits kinking, they easily maneuver through narrow spaces and curves, providing surgeons with precision, reliability, and confidence in every procedure.

Zebcor[™]
Zebra Guidewire

Ribrid[™]
Hybrid Guidewire

Hydroease[™]
Hydrophilic Guidewire

Meriskate[™]
PTFE Guidewire





ZebcorTM

Zebra Guidewire

Features


- Excellent surface properties makes it super lubricious.
- It is simpler to see the movement of the guidewire because of the distinctive white and blue stripes.
- The super flexible nitinol core's exceptional torsional stiffness prevents kinking.
- The tapered tip shape allows the guide wire to navigate narrow spaces and curves more easily.

Size Matrix

STRAIGHT	J
0.035"inch x Length: 150CM	0.035" inch x Length:150CM
0.028" inch x Length:150 CM	

Use

It facilitate the placement of devices during urological procedures by providing a smooth, flexible pathway.



Ribrid™

Hybrid Guidewire

Features


- Combines a stiff shaft for pushability with a flexible tip for safe navigation.
- Distal hydrophilic coating ensures superior lubricity in wet conditions.
- Non-stick shaft coating enhances grip, torque control, and insertion speed.
- Highly fluoroscopic-visible for precise positioning and monitoring.

Size Matrix

STRAIGHT	J
0.035" Inch X Length:150CM	0.035" Inch X Length:150CM

Use

The hybrid guidewire aids navigation and access in urological procedures.



HydroeaseTM

Hydrophilic Guidewire

Features


- Activates in wet conditions for ultra-low friction, enabling smooth navigation through strictures and complex anatomy.
- One-Pass Access with Smooth Trackability.
- Kink-Resistant and Memory-Stable.
- Hydrophilic layer stays active throughout the procedure for continuous smooth glide.

Size Matrix

STRAIGHT	CURVED
0.032"inch x Length: 150 CM	0.032"inch x Length: 150CM 0.035" inch x Length: 150CM
0.035"inch x Length: 150 CM	
0.038"inch x Length: 150 CM	

Use

The hydrophilic guidewire aids smooth navigation in urinary procedures.



MeriskateTM

PTFE Guidewire

Features

- Hydrophilic coating minimizes friction to the greatest extent feasible.
- The flexible tip provides atraumatic advancement.
- Excellent torquability is intended for precise control and placement.

Size (Fr)

0.032", 0.035"

Coating


Hydrophilic Coating

Working Length

150cm

Tip

1. Straight Tip
2. Curved Tip



The Next Era of **Urology**



Dilator

Gentle Expansion, Effortless Precision

Like the artful rhythm of a gymnast's performance, where every movement is deliberate and fluid, our dilators are crafted for smooth and controlled expansion. With optimal flexibility and precise dilation, they navigate anatomical pathways effortlessly, minimizing trauma while maximizing procedural efficiency—empowering surgeons with confidence at every step.

Dilato[™]
Fascial Dilator Set

Dilatatrix[™]
Urethral Dilator

EasyDi[™]
Ureteral Dilator



DilatoTM

Fascial Dilator Set

Features


- Facilitate smooth surface for atraumatic entry.
- High durability.
- Effective stricture dilation improves outcomes with fewer complications ad greater patient comfort.

Size Matrix

L:18 CM	L:22 CM
Dilator: 8-18 Fr	(6 to18 FR)
	(6 to16 FR)

Use

It aids passage of nephrostomy or Malecot catheters.



DilatrixTM

Urethral Dilator

Features



- Maintains shape while adapting to anatomical curves.
- Tapered end with outstanding durability and resilience.
- Biocompatible material results in better experience.

Size Matrix

- | | | |
|------------------------|----------------------|----------------------|
| • 8Fr x Length : 37CM | • 14Fr X Length:37CM | • 20Fr x Length:37CM |
| • 10Fr x Length : 37CM | • 16FrX Length:37CM | • 22Fr x Length:37CM |
| • 12Fr x Length : 37CM | • 18FrX Length:37CM | • 24Fr x Length:37CM |

Use

It dilates the urethra before ureteroscopy or stone removal.



EasyDiTM

Ureteral Dilator

Features




- Smooth dilator surface ensures low-friction insertion with minimal patient trauma.
- Tapered end with outstanding durability and resilience.
- Biocompatible material results in better experience.

Size Matrix

- 6FR X Length:60CM
- 8FR X Length:60CM
- 10FR X Length:60CM
- 12FR X Length:60CM
- 14FR X Length:60CM
- 16FR X Length:60CM

Use

The Ureteral Dilator widens the ureter before ureteroscopy or stone removal.



Scope

Clarity in Vision, Precision in Action

Just as binoculars bring distant details into sharp focus, our ureteroscopes provide exceptional imaging for precise diagnostics and treatment. Equipped with a high-performance CMOS camera and advanced LED illumination, they deliver high-definition visualization, ensuring accurate navigation, enhanced clarity, and superior procedural outcomes—empowering surgeons with confidence at every step.

ri-Uniview[™]
Flexible Fiberoptic Ureteroscope

riUroSnap[™]
Image Processor without Screen



ri-UniviewTM

Flexible Fibreoptic Ureteroscope

Features

- The digital flexible ureteroscope has both the apex and the insertion tube maintained at 7.5Fr & 6.3Fr, with the apex capable of bending in four directions (275 degrees up and down, 90 degrees left and right) for imaging.
- The tip, with a 3.5Fr inner diameter, features a CMOS camera and an LED.

Size Matrix

L:65 CM	L:68 CM
8.5 Fr	6.3 FR
	7.5 FR

Use

The Ureteroscope System visualizes and treats urinary tract areas via transurethral or percutaneous access.



riUroSnapTM

Image Processor without Screen

Features

- Delivers ultra-clear, high-contrast visualization with true color fidelity.
- Equipped with multiple video output ports for connecting additional screen.
- Enables data recording with USB drive storage capability

Video Resolution

1920 x 1080 P

160k



The Next Era of **Urology**

Laser and Special Category Performance in Motion, Precision in Action

Just as binoculars bring distant details into sharp focus, our ureteroscopes provide exceptional imaging for precise diagnostics and treatment. Equipped with a high-performance CMOS camera and advanced LED illumination, they deliver high-definition visualization, ensuring accurate navigation, enhanced clarity, and superior procedural outcomes—empowering surgeons with confidence at every step

ProNeph™

Percutaneous Nephrolithotomy (PCNL) Package

MeriFRAG™+

Holmium Laser Portable System



MeriFRAG+

Holmium Laser Portable System

MeriFRAG 40 WATTS HOLMIUM:LASER FOR LITHOTRIPSY

- MeriFRAG surgical laser is a Holmium:
- YAG (HO:YAG) laser that is designed for Lithotripsy, Coagulation and tissue Cutting/ Ablation. The Laser with a wavelength of
- 2100nm is highly absorbed by water and biological tissue provide excellent cutting, ablation and hemostatic properties.
- The surrounding tissue sustains minimal damage as a result of the restricted radiation penetration (0.3–0.4 mm).

MeriFRAG HOLMIUM LASER SYSTEM

- Effective Lithotripsy.
- Minimal Penetration
- Soft Tissue Surgery.
- Compact Design and Light Weight.
- All-Round Solution for all the Stones.

Application Details

MeriFRAG can be used to execute a variety of medical specialties' procedures, including soft tissue's incision, excision, resection, ablation, vaporization, coagulation and hemostasis as well as stone lithotripsy

Different Modes of Lithotripsy

1. Fragmentation Mode.
2. Popcorn Mode.
3. Dusting Mode.

Technical Specifications

Output Power (Max)	40 WATTS
Pulse Frequency	1 – 30Hz
Pulse Energy	0.1 – 6.0J
Laser Type	Holmium Laser (Ho:YAG)
Wavelength	2.1 μ m
Pulse Duration	50 - 1500 μ s
Aiming Beam	532nm < 1mW, can be fine-tuned
Cooling system	Internal Cooling System
Fibers	Reusable 200, 276, 365, 550, 800, 1000 μ m
Main Supply	220V+22V/ 50Hz + 1Hz Available
Dimension and Weight	49 cm(H) x 46 cm(W) x 56 cm(L) – 45kg



ProNephTM

Percutaneous Nephrolithotomy (PCNL) Package

Features

To safely create a percutaneous tract into the kidney, dilate the access, and facilitate stone removal using appropriate instruments.

Package Includes:

- 1. Fascial Dilator Set
- 2. PTFE/Hydrophilic Guide wire
- 3. Puncture Needle
- 4. Peel Away Sheath
- 5. Nephrostomy Tube

Sizes

Peel away sheath:

L:19CM	12Fr	14Fr	16Fr	18Fr
--------	------	------	------	------

Nephrostomy tube :

L:18CM	12 FR	14 FR	16 FR	18 FR	20 FR	22 FR	24 FR
L:22CM	12 FR	14 FR	16 FR	18 FR	20 FR	22 FR	24 FR
L:22CM	12 FR	14 FR	16 FR	18 FR	20 FR	22 FR	24 FR

Puncture needle:18 Goval tip Length: 20 CM

Hydrophilic/PTFE Guidewire : 0.035"inch Length:150CM


Fascial Dilator Set :

L:18CM	8 FR	10 FR	12 FR	14 FR	16 FR	18 FR	22 FR
L:22CM	8 FR	10 FR	12 FR	14 FR	16 FR	18 FR	22 FR

Inner Diamete (ID) (mm): 4.30, Outer Diameter (OD) (mm): 5.75

Use

The PCNL package used for kidney stone removal PCNL Procedure.



Leading the Way in Advanced Urological Solution

Urology Solutions Upcoming Portfolio

- Urine Bag
- Thulium Laser Machine
- Electro-hydraulic Lithotripter
- Urodynamic System
- Uroflowmeter Meter
- Artificial Urinary Sphincter
- Artificial Penile Prosthesis
- Triple J Stent
- Chiba Needle
- CIC Catheter
- Re-entry Malecot Catheter
- Bladder Evacuator
- Irrigation Bulb
- Meatal Dilator Male
- Meatal Dilator Female
- Semi-auto Biopsy Needle
- Automatic Biopsy Needle
- Malecot Catheter
- Tur Set
- Urine Meter
- Irrigation Syringe
- Silicon Foley Catheter
- Nephroscope
- Cystoscope
- Resestoscope
- Endopyelotomy Stent
- Cytoscopic Injection Needle
- Silicon Double Pigtail Ureteral Stent
- Collins Knife
- Cutting Loops
- Cold Knife
- Bubbee Electrode
- Bi Polar Button
- Bipolar Cutting Loops
- Roller Electrodes
- Cath Tap
- Urinary Foley Catheter with Temperature Sensor
- Introducer Needle
- Grasping Forceps Bi Prong Rigid
- Grasping Forceps Bi Prong Flexible
- Grasping Forceps Tri Prong
- Nephrostomy Balloon Dilator Catheter
- PCNL Step Dilator
- Screw Dilator

- Urethral Self Dilator
- Single Step Ureteral Dilator
- Percutaneous Nephrostomy Set
- Amplatz Renal Dilator Set
- Delta Basket
- Calicone
- Helical Basket
- Nephrostomy Balloon Dilator
- Injection Needles with Side Holes
- Reelief CI Catheter (Male)
- Reelief CI Catheter (Female)
- Alligator Peanut Forcep
- Ureteroscopy Irrigation System
- Disposable Two Parts Trocar Needle
- Universa Loop Nephrostomy Perc Drainage Cathter
- Malecot Supra Perc Drainage Catheter
- High Pressure Balloon Cathter
- Single Action Pumping System
- Water Vapor Therapy
- Cystoscope Sheath
- Flexible Cystoscope
- Stent Removing Forceps
- Stone Punch SET
- PCN Catheter SET
- Prostatic Stent
- AQUA – SUI Solution Kit
- Dual Lumen Catheter
- Rectal Balloon Catheter
- Urethral Sling
- Path Finder
- Penile Clamp
- Morcescope Set
- Prostate Enucleation Solutions
- Urolift
- Prostatic Stents
- Metal Stent
- Flexible Nephroscope
- MINI PCNL Nephroscope
- Closed end Ureteral Catheter
- Urethral Stent
- Rapid Release-Ureteral Access Sheat

Meril

More to Life

Manufactured by

Meril Medical Innovations Pvt. Ltd.

Meril Park 2, F5, F6 Survey No. 879, Muktanand Marg, Chala,
Vapi, Valsad, Gujarat, 396191, India

Tel: +91 260 240 8000 Fax: +91 (260) 240 8025

E-mail: askinfo@merillife.com Website:

<http://www.merillife.com>

